

December 10, 2021

AUTOCLAVES

OPTION 1

**MANUAL AUTOCLAVES**

**R26000.00 excluding VAT**

**PORTABLE SA-232 SERIES**



**CODE : ASA232**

**FEATURES :**

- Inexpensive, reliable, safe, easy operation and maintenance
- High efficiency for sterilizing wrapped/unwrapped/liquid/utensils, etc.
- Stainless steel sheathed heater for easy maintenance
- Versatile programs for selection:
  - Sterilization pressure adjustable from 0.9 – 2.1 kg/cm<sup>2</sup> assoc. with 115°C – 134°C
  - Preset sterilization time at 4 / 15 / 30 minutes depending on sterilization substances
- Microprocessor control operates sterilization automatically
- Sterilization starts timing only once set temperature is reached
- Complies with quality certifications and specifications of 93/42/EEC

**SPECIFICATIONS :**

<b>MODEL :</b>	<b>SA-232 Portable</b>
<b>Overall Dimensions</b>	330mm x 420mm x 400mm
<b>Capacity</b>	16 ℓ
<b>Temperature / Pressure</b>	Unwrapped / wrapped / liquid : 122°C ; 1.4kg per cm <sup>2</sup>
<b>Control System</b>	Pressure control
<b>Power Consumption</b>	220-240V: 0.9kW; 4Amp.
<b>Auto Drying</b>	N/A
<b>Pre-vacuum, Dry-vacuum</b>	N/A
<b>Indicators</b>	Yes
<b>Temperature Display</b>	Temperature gauge
<b>Sterilization Temp / Time</b>	Wrapped: 121°C - 30 minutes

	Unwrapped: 121°C - 15 minutes
<b>Standard Accessories</b>	<ul style="list-style-type: none"> <li>✓ Sterilization box x 1</li> <li>✓ Tray x 1</li> </ul>
<b>Optional Accessories</b>	<ul style="list-style-type: none"> <li>✓ Tray frame x 1</li> <li>✓ Tray x 3</li> <li>✓ Tray holder x 1</li> <li>✓ Spring x 2</li> </ul>
<b>Safety Devices</b>	<ul style="list-style-type: none"> <li>✓ Over-heat shutout</li> <li>✓ Safety valve</li> <li>✓ Release valve</li> <li>✓ Emergency exhaust valve</li> <li>✓ Door safety holder</li> <li>✓ Pressure control switch</li> </ul>
<b>Quality Standard</b>	CE-MDD

---

OPTION 2

## MANUAL TABLE TOP AUTOCLAVE SA-232X SERIES

**R39500.00 excluding VAT**



**CODE : ASA232X**

**FEATURES:**

- Inexpensive, reliable, safe, easy operation and maintenance
- High efficiency for sterilizing wrapped/unwrapped/liquid/utensils, etc.
- Stainless steel sheathed heater for easy maintenance
- Versatile programs for selection:
  - Sterilization pressure adjustable from 0.9 – 2.1 kg per cm<sup>2</sup> assoc. with 115°C – 134°C
  - Preset sterilization time at 4 / 15 / 30 minutes depending on sterilization substances
- Microprocessor control operates sterilization automatically
- Sterilization starts timing only once set temperature is reached
- Complies with quality certifications and specifications of 93/42/EEC

**SPECIFICATIONS :**

MODEL :	232X Table top
Actual Dimensions	500mm x 400mm x 530mm
Capacity	16 ℓ
Temperature / Pressure	Unwrapped / wrapped / liquid : 115°C -134°C ; 0.9-1.4kg per cm <sup>2</sup>
Control System	Pressure control
Power Consumption	220-240V : 1.5kW ; 6.5Amp
Auto Drying	N/A
Pre-vacuum, Dry-vacuum	N/A
Indicators	Yes
Temperature Display	Temperature gauge
Sterilization Temp / Time	Wrapped: 121°C - 30 min; 134°C - 15 min Unwrapped: 121°C - 15 min; 134°C - 4 min
Standard Accessories	✓ Sterilization box x 1 ✓ Tray x 1
Optional Accessories	✓ Tray frame x 1 ✓ Tray x 3 ✓ Tray holder x 1 ✓ Spring x 2
Safety Devices	✓ Over-heat shutout ✓ Safety valve ✓ Release valve ✓ Emergency exhaust valve

---

**Ecotao Enterprises cc (VAT no: 4590247682)**

39 Olienhout Avenue, Platteklouf3, Cape Town, 7500, South Africa

| Tel: 27-21-9115835 | Skype: ecotao | email: aqua@ecotao.co.za | Mobile: 27-823903989|

---

	✓ Door safety holder
	✓ Pressure control switch
<b>Packaging Dimensions</b>	640mm x 590mm x 480mm
<b>Packaging Weight</b>	40kg
<b>Quality Standard</b>	CE-MDD